Laboratory Techniques Overview

Na	nme:	Section:	Partners:	
cla we tec	ass each week. Each technique is refere eek you will be required to conduct seve	nced throughout the ral activities to far ory. Please read thr	nsible for reading and understanding before current weeks laboratory experiment philiarize yourself with the five most comough and complete each activity. Once	procedure. This nmon
I.	A. One common technique for fin water in a graduated cylinder is final volume of the water in the final volume will give the volu Using the volume data listed be just like it is shown in the Tech	ons in your report" Fork just like it is ding the volume of segraduated cylind are of water displayed on, calculate the	and "Significant figures" techniques. Poshown in the Techniques Document. If an object is by displacement. The initial piect is then placed into the graduated cyper is recorded. Subtracting the initial voiced by the object. This is the volume of volume of the object for the two trials.	al volume of linder and the lume from the f the object.
	Table 1: Volume Data	Twist 1	Trial 2	
	Initial Volume	Trial 1 51.0 mL	50.0 mL	
	Final Volume	62.5 mL	60.5 mL	
	Formula:			
	Work:			
	Answer:			
	B. Another common calculation is Formula:	n lab is to find an a	verage. Find the average of the two vol	umes from IA.
	Work:			

Answer: _____

C. The volume of a rectangular block can be found by calculation. Using the given data, find the volume of the rectangular block.

Table 2: Rectangular Block Data

Length	Width	Height	
12.25 cm	7.90 cm	3.68 cm	
Formula:			
Work:			

Answer:	

II. Reading a Balance: (10 points)

Read through the "Using and analytical balance" technique. Find the station with two objects that you need to weigh next to a balance. Weigh each object and correctly record the weight on your lab report sheet.

- A. Weigh the aluminum rectangular prism: ______
- B. Weighing out solid chemicals requires a few extra steps. One common technique is to place a weigh boat on the balance and then tare the weigh boat. Obtain a weigh boat from next to the balance and press the tare button. Make sure all of the glass doors are closed tight. Using a scoopula, weigh out approximately 1-g of NaCl. Try to get as close to 1-g as you can without going over. If you realized you scooped too much, the correct procedure is to place the excess into a waste beaker, NEVER back into the original container. Placing excess back into the original container can introduce impurities and contaminate the entire container. This is true for liquids as well.

Record the exact weight of the NaCl from the balance here:

Empty your weigh boat into the waste beaker next to the balance and clean the balance area. There is a brush next to the balance to brush off any solid that spilled on or around the metal tray. This is important to do each and every time you use the balance since some chemicals can react with the materials the balance is composed of and over time this does destroy our balances.

III. Reading a Meniscus: (10 points)

Read through the "Reading a meniscus in various pieces of glassware" technique. Find the station set-up with common pieces of glassware and record the volume of the meniscus for each below:

		First Attempt:	Second Attempt: (only if required)
A.	100-mL graduated cylinder		
B.	10-mL graduated cylinder		
C.	25-mL buret		

Check your answers with your instructor. You may need to go back and try this again.

- D. Go to the hood and use the pump to dispense 50-mL of water into a 150-mL beaker. The pump is pre-set to 50-mL so do not adjust the setting. Look at the increment marks on the pump. What decimal place should the pump volume be recorded to?
 - a. ones place
 - b. tenths place

Explain your findings:

- c. hundredths place
- d. thousandths place

Read the volume on the 150-mL beaker and record your answer here to the correct significant figures			
Does the volume on the beaker match the volume on the pump? Yes or No.			

IV. Operating a Bunsen Burner: (10 points)

Find where the Bunsen burner and striker are stored for your lab bench. Pull it out and check the hose for cracks. Follow the techniques document "Using a Bunsen burner" to create a 2-inch flame with an inner and outer blue cone. Once each student at your desk has successfully operated the Bunsen burner let your instructor know you are ready for his or her signature. You will be required to demonstrate your technique to your instructor before obtaining their signature. Always turn off the Bunsen burner at the gas spigot on the bench, never by rotating the dial or shaft. Turning the dial or shaft does not turn off the gas!

V. Handling Chemicals: (10 points)

Read through the "Handling chemicals" section of the techniques document. Understanding how to handle your waste is essential for a general chemistry laboratory. All solid and liquid waste must be placed in an appropriately labeled waste container, never down the sink. The waste containers for this course will always be located in the hood. Sometimes two different types of waste containers may be necessary depending on the experiment and the chemical waste generated.

- A. Locate the two different waste containers in the hood.
- B. Draw the waste containers in the space below. Include the waste container label on your drawing with the exact contents of the waste containers.